The **SMART BAG®** Disposable Bag Valve Mask Resuscitator is a single use device for manual ventilation. The **SMART BAG®** offers accurate controlled ventilation that significantly reduces the risk of gastric insufflation. The patient's airway pressure is maintained at an optimally low level due to the **SMART®** technology in the bag that responds to both the rescuer's squeeze of the bag and the patient's own respiratory condition. It is designed for use in various clinical settings to provide respiratory support to the respiratory distressed or non-breathing patient.

The Child resuscitator comes standard with a Pressure Relief Valve set at 40 cm H₂O which can be overridden by the operator. A pressure relief system, set at either 40 or 60 cm H₂O, is available for the Adult resuscitator. When used with a supplemental oxygen supply, the device may provide 100% oxygen on every squeeze of the bag.

PREPARATION FOR USE

[1]Inspect the SMART BAG® resuscitator to ensure that all components are present, properly assembled and fully functional.

[2] Test for leaks by occluding the patient port completely squeezing the bag. Any leaks in the system may prevent the delivery of a sufficient volume to the patient.

[3] Squeeze and release the **SMART BAG®** hard a few times to ensure that air is moving through the valve system to the mask. The **SMART®** Valve in the neck of the bag should move freely indicating increased airway pressure and you should notice an immediate increase in bag tension (stiffness).

[4] Gently squeeze and release the **SMART BAG®** a few times to ensure that the bag tension is reduced and the **SMART®** Valve in the neck or the bag does not move forward when you gently squeeze. This provides confirmation that the airway pressure will be kept to the minimum required for adequate ventilation to occur while reducing the risk of gastric insufflation.

[5] If using supplemental oxygen, attach the reservoir system to the bag refill port and ensure that the oxygen tubing is attached to an oxygen source with a flow rate of at least 15 L/min. Ensure that the collapsible reservoir system is fully extended to allow maximum oxygen storage.

DIRECTIONS FOR USE

[1] Select the appropriate **SMART BAG®** resuscitator model for the size of patient to be ventilated.

[2] Ensure that the patient's airway is clear of any obstructions and remains open by properly positioning the patient's head in accordance with local protocols.

[3] Maintain a proper mask-to-face seal with one hand by lifting the chin upward with the last three fingers of the hand. Keep the index finger and thumb on top of the mask to form a tight seal around the patient's mouth and nose.

[4] Gently squeeze the **SMART BAG®** with the other hand until the chest rises, then release. Ventilate the patient with a steady squeeze and release of the **SMART BAG®** sufficient time between ventilations to allow for full emptying of the patient's lungs.

[5] If the child **SMART BAG®** is being used and the Pressure Relief override is required to be applied, depress the Pressure Relief Button and rotate 90° to lock in place. To unlock simply rotate the button until the arrow lines up with the arrow on the patient valve and release.

[6] Safely dispose of the **SMART BAG®** after use in accordance with local protocols.

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