



O-Two Mini Ventilation Training Analyzer

Manual Ventilation Training

Teaches the provider to eliminate
inadvertent hyperventilation

O-Two Mini Ventilation Training Analyzer

Healthcare providers being trained in bag-valve-mask ventilation are generally taught to ventilate until the chest rises. However, studies have shown that the delivery of incorrect rates and volumes is common.

Inadvertent hyperventilation (the accidental delivery of an excessive minute volume) has been shown to cause increased cerebral ischemia, decreased coronary perfusion pressure and gastric insufflation. In the unprotected airway, if the lungs are over-inflated or the delivered flow rate is too high, the high airway pressure generated will cause air to enter the stomach. In real life this can only be seen when the patient's abdomen distends and/or they begin to vomit.

The mini ventilation training analyzer has been developed to demonstrate the correct ventilation technique when using a bag-valve-mask resuscitator (BVM) and the complications of gastric insufflation associated with incorrect technique.

Based on a "normal" adult lung capacity, compliance, airway resistance and lower esophageal

sphincter pressure (LESP) of 19 cmH₂O, the mini ventilation training analyzer responds to the Healthcare provider's squeeze of the BVM like a real patient.

If Healthcare providers squeeze the BVM too quickly or forcefully, the high pressure they generate (due to the high air flow from the BVM) will cause air to enter the "stomach" through the calibrated lower esophageal sphincter.

Designed for manual ventilation training at all levels of rescuer — doctors, nurses, respiratory therapists, paramedics, EMTs, fire fighters, first responders, CPR instructors and students, the mini ventilation training analyzer is a unique training tool that teaches real-life skills while identifying real life problems.

It can also be used for training in the proper use of manually triggered, gas-powered resuscitators.



Carry case for 01TA1200-MO

Features

- Small, lightweight and easily portable
- Simple to assemble
- Calibrated "lower esophageal sphincter"
- Simulated stomach
- Calibrated test lung with compliance band and normal airway restriction
- 22/22 mm BVM connector
- 30 mm "stomach" connector
- Unique lung connector
- Instructions for use printed on top of the training analyzer

Ordering Information

01TA1000 Mini Ventilation Training Analyzer with Calibrated Test Lung, Stomach and Lower Esophageal Sphincter Valve

01TA1200-MO Mini Ventilation Training Analyzer
Complete training kit in softpack carrying case. Includes SMART BAG MO controlled ventilation BVM and Easy Grip® standard BVM.



7575 Kimbel Street, Mississauga, ON, Canada, L5S 1C8
Telephone +1 905 677 9410 N.A. Toll Free +1 800 387 3405
Facsimile +1 905 677 2035 Web otwo.com
Email resuscitation@otwo.com

SMART BAG and Easy Grip are registered trademarks of O-Two Medical Technologies Inc.