

O-Two Self-Study Guide

**e500 Transport Ventilator
Ventilation Modes**



E500 VENTILATION MODES


The e500 ventilator has 5 ventilation modes:

Manual Ventilation, Synchronized Intermittent Mandatory Ventilation (SIMV), CPAP, Mask CPR, intubated CPR. In addition the ventilator provides a “Quick Start” mode when the ventilator is turned on.

1. Quick Start Mode

The default ventilation mode for “Quick Start” is SIMV.

There are 12 tidal volume/rate options on the startup screen:




Rotate the control knob to move the cursor over the required patient size select and confirm by pressing the control knob.

Note: If a patient size selection is not made and confirmed within 20 Seconds of start up, the e500 Ventilator will automatically commence ventilation in the 500/10 setting. Following startup it is easy to make setting changes to ventilation parameters, alarms or to select other ventilation modes.

The pre-set “Quick Start” defaults are based on the American Heart Association Resuscitation Guidelines recommended Tidal Volume/Ventilation Frequencies for the patient sizes indicated.

2. Synchronized Intermittent Mandatory Ventilation - (SIMV)


The SIMV mode allows the patient to breath spontaneously. If assistance is needed the vent will assist. If NO patient effort during trigger window - Ventilator will initiate Mandatory Ventilation.



Screen display for SIMV mode

Parameter	Range	Default
Tidal Volume/ Rate	100/20,150/20, 200/20, 250/20, 300/10, 350/10, 400/10, 450/10, 500/10, 550/10, 600/10, 650/10	500/10 or selected by user during start-up
I:E ratio	1:2	1:2
PEEP	0, 4, 5, 6, 7, 8, 9, 10, 12, 15, 18, 20	5
Trig.	3 L/min	3 L/min
O ₂ %	(100% or 60%)	100%
P max.	10, 15, 20, 25, 30, 35, 40, 45, 50, 55, 60, 80	30
P min.	0, 2, 3, 4, 5, 6, 7, 8, 10, 12, 15, 20	3
Manual	Refer to Manual and I-Hold section	ready


Range and Ventilation Settings for SIMV.



SIMV Waveform

3. Constant Positive Airway Pressure (CPAP)

In CPAP mode, the ventilator will deliver a continuous flow rate to generate a continuous positive pressure in the airway. The ventilator uses the control valve to maintain the CPAP level.




Screen display for CPAP mode

Default Ventilation Setting - CPAP

Parameter	Range	Default
CPAP (cm H ₂ O)	4, 5, 6, 7, 8, 9, 10, 12, 15, 18, 20	Carry over from PEEP of previous mode
P max.	10, 15, 20, 25, 30, 35, 40, 45, 50, 55, 60, 80	Carry over from previous mode
P min.	0, 2, 3, 4, 5, 6, 7, 8, 10, 12, 15, 20	Carry over from previous mode
O ₂ %	(100% or 60%)	Carry over from previous mode
T APNEA (sec)	10, 15, 20, 25, 30, 40, 45, 50, 55, 60	20
Vt/f (A) (ml/BPM)	100/20, 150/20, 200/20, 250/20, 300/10, 350/10, 400/10, 450/10, 500/10, 600/10, 650/10	Carry over from previous mode

Note: CPAP is pre-set with 20 second APNEA backup - settings for Rate BPM, Tidal Volume (Vt), I:E ratio, and Minute Volume (Mv) location is set at the bottom of the screen. The default settings can be adjusted by the operator.



CPAP waveform

4. CPR Modes

The CPR Mode consists of timed chest compression audible prompts coupled with automatically delivered breaths for both intubated and mask ventilated patients. The screen displays a visual animated display to guide the health care provider while performing CPR. Oxygen is delivered automatically at 100%.


Default Ventilation Setting- CPR

Parameter	Range	Default
Tidal Volume	100,150, 200, 250, 300, 350,400, 450, 500, 550, 600, 650	Carry over from previous mode
P max.	10, 15, 20, 25, 30, 35, 40, 45, 50, 55, 60, 80	Carry over from previous mode
P min.	0, 2, 3, 4, 5, 6, 7, 8, 10, 12, 15, 20	Carry over from previous mode
O ₂ %	100%	100%




Mask Ventilated Patient Basic Airway - 30:2

**With each chest compression the hands are highlighted;
the mask highlights in blue when a breath is delivered.**




Mask CPR Ventilation Waveforms



**Intubated Airway Ventilation
Advanced Airway - Asynchronous**

With each chest compression the hands are highlighted; the ET tube highlights in blue when a breath is delivered



Intubated CPR Ventilation Waveforms

O-Two Medical Technologies Inc.

45A, Armthorpe Road, Brampton, Ontario, Canada L6T 5M4

Toll Free: **(800) 387 3405** Tel: **(905) 792 OTWO (6896)** Fax: **(905) 799 1339**

Email: resuscitation@o-two.com Website: www.o-two.com

NOTE: Refer to the product manual for further definitions and descriptions of function.

Rev 0. 23_08_2017