

# CPAP — VS — BPAP

Quick Guide



## Step 1 Check if the patient is suitable for (CPAP or BPAP)

- ✓ Awake, alert, and cooperative
- ✓ Spontaneously breathing, able to protect airway
- ✓ Hemodynamically stable (no severe hypotension)
- ✗ No facial trauma, vomiting, or massive secretions


## Step 2 Identify the main problem (oxygenation vs ventilation)

### ● OXYGENATION ONLY

- SpO<sub>2</sub> (<92%) despite O<sub>2</sub> therapy.
- ETCO<sub>2</sub> is normal or low.
- Clinical signs: cyanosis, crackles, pulmonary edema, pneumonia, near drowning.

→ Use **CPAP**

- Start **CPAP** level at 5-10 cmH<sub>2</sub>O (range 4-20).
- Set **O<sub>2</sub>** at 60% or 100%, titrate to target SpO<sub>2</sub>.
- Set Apnea backup **T-APNEA** (Default: 20s). If Apnea occurs, the ventilator switches to A/CV mode.
- Monitor comfort, chest rise, work of breathing, and oxygenation.




### ● OXYGENATION + VENTILATION

- Low SpO<sub>2</sub> and rising ETCO<sub>2</sub> >45 ↑
- Signs of fatigue (↑ RR, shallow breaths, accessory muscle use).
- Example: Acute exacerbation of COPD, Severe Asthma, and Cardiogenic pulmonary edema

→ Use **BPAP** (Biphasic Positive Airway Pressure)  
CPAP + PSV

- **EPAP**: Set the **CPAP** baseline (Ex: 5 cmH<sub>2</sub>O).
- Add **PSV** 5-15 cmH<sub>2</sub>O (titrate by 2-3).
- **IPAP** is determined by summing the **CPAP** and **PSV** values.
- Keep Peak Airway Pressure **PAWpeak** ≤30 cmH<sub>2</sub>O (set **Pmax** alarm >30 cmH<sub>2</sub>O).
- (**PAW peak** in this mode will be the sum of **CPAP+PSV** values)



## Step 3 Escalate or switch if the patient worsens

- If on **CPAP** → develops ↑ ETCO<sub>2</sub> / ↑ work of breathing/ fatigue  
→ **Switch to BPAP**
- If on **BPAP** → deteriorates or loses consciousness  
→ **Prepare for intubation**

### DON'T DELAY AIRWAY IF REQUIRED

Treatment Failure = worsening hypoxemia/hypercapnia, rising ETCO<sub>2</sub> despite support, decreasing LOC, inability to tolerate mask, refractory hypotension, or absent improvement within ~10-20 min → move to advanced airway per protocol.

## Capnography Quick Guide

Use **Capnography (ETCO<sub>2</sub>) with Pulse Oximetry (SpO<sub>2</sub>)** to guide mode selection.

This applies only to awake, spontaneously breathing patients appropriate for non-invasive ventilation (**CPAP** or **BPAP**).

ETCO <sub>2</sub> (mmHg)	Interpretation	Oxygenation Status (SpO <sub>2</sub> )	EMS Action
<b>35-45</b>	Normal ventilation	Low SpO <sub>2</sub> only	Use <b>CPAP</b> – recruit alveoli, improve oxygenation
<b>&gt;45</b> (hypercapnia)	CO <sub>2</sub> retention / impaired ventilation	Low SpO <sub>2</sub> or normal	Use <b>BPAP</b> – add pressure support (PSV) to assist ventilation
<b>&gt;50</b>	Significant CO <sub>2</sub> retention	Low SpO <sub>2</sub> or normal	Start urgent <b>BPAP</b> ; if no improvement, prepare for intubation
<b>Rising trend</b> (+5 - 10 mmHg in minutes)	Worsening ventilation	Any SpO <sub>2</sub> status	Escalate ventilatory support (invasive ventilation), monitor closely

## Step 4 Monitor continuously

- SpO<sub>2</sub>, ETCO<sub>2</sub>, RR, HR, BP
- Check mask seal and comfort every few minutes
- Watch for hypotension, gastric distention, or decreased Level of Consciousness (LOC).

If ETCO<sub>2</sub> >45 mmHg AND patient shows fatigue or altered LOC, treat as hypercapnic respiratory failure

→ Use **BPAP** (if no contraindications.)

## Clinical Signs Supporting Hypercapnia

- Increased work of breathing (accessory muscles, retractions)
- Decreasing respiratory rate with shallow breaths
- Altered mental status or drowsiness (CO<sub>2</sub> narcosis)
- Fatigue, inability to speak full sentences

Disclaimer: This document is for informational and educational purposes only and is not a clinical guideline, protocol, policy, or standard of care. It does not replace clinical judgment or local/regional EMS directives.

[www.otwo.com](http://www.otwo.com)

45A Armthorpe Road, Brampton, ON, Canada, L6T 5M4  
Telephone: +1 905 792-OTWO (6896) N.A. Toll Free: +1 800 387 3405  
Facsimile: +1 905 799 1339 Email: [resuscitation@otwo.com](mailto:resuscitation@otwo.com)

**o\_two** controlled™  
ventilation