

Manual vs. Mechanical Ventilation in Prehospital Setting

Alba Marin, MD | O-Two Clinical Specialist

Summary

Even though bag-valve-mask system is one of the most common devices used to provide ventilation during transports, it has shown that is complicated, has variable results, makes it difficult to reach the standards in terms of tidal volume and respiratory rate, and may expose to overpressure and thoracic overinflation. Due to the variability in respiratory parameters and changes in hemodynamics during manual ventilation, mechanical ventilatory support with portable devices have become the preferred method for transporting patients in prehospital environments. Mechanical ventilators offer a safe alternative to manual ventilation and allow healthcare providers to deliver consistent care and controlled tidal volumes at a determined rate. Patients can benefit from early initiation of non-invasive treatment in prehospital setting to avoid intubation and to improve patient outcomes.

Airway management and ventilation assistance are some of the most important aspects when managing out-of-hospital medical emergencies^{1,2}. They directly affect patients' potential for recovery, promote oxygenation, and may protect against aspiration depending on the approach. Optimal controlled ventilation is crucial to optimize inspiratory time, expiratory time, and airflow³.

Numerous techniques and devices are available to deliver oxygen-enriched air to patients during resuscitation⁴ and acute respiratory failure. Historically endotracheal intubation has been considered the gold standard for airway management. Nowadays, the primary objective in prehospital environments is to assure ventilation of the patient, invasively or non-invasively, during transport to an emergency department. Options for airway management and ventilatory support involve different levels of invasiveness and complexity that require different technologies and expertise³.

Manual ventilation during transport

The bag-valve-mask (BVM) is one of the most common devices used to provide ventilation, although the American Heart Association (AHA)

ranks BVM devices lower in preference than other ventilation methods, such as emergency and transport ventilators³⁻⁵. **BVM resuscitation makes it difficult to reach the standards** in terms of tidal volume (Vt) and respiratory rate (RR)^{2,6,7} and may expose to overpressure and thoracic overinflation^{8,9}. Ventilation rates during the on-field application of cardiopulmonary resuscitation (CPR) by well-trained emergency medical services (EMS) personnel are far in excess of those recommended¹⁰.

The **AHA Guidelines for Cardiopulmonary Resuscitation and Emergency Cardiovascular Care** recommend ventilation rates of eight to ten breaths per minute or two ventilations every 30 compressions, and tidal volumes between 500-600 mL¹¹. The recommendation is to deliver each **rescue ventilation over 1 s**, give a sufficient tidal volume (Vt) **to produce a visible chest rise**^{6,11,12}. This given volume of ventilation is adequate for oxygenation and minimizes the risk of gastric inflation¹². However, both the provider and observer cannot be certain of the amount of Vt delivered, because it varies depending on the operator and squeeze method. The Vt delivered by a standard manual resuscitator shows large variations despite being performed by well-trained EMS providers^{3,10}.

Moreover, **manual ventilation is not recommended to be performed by a single provider**^{12,13}. It requires two hands: one hand to secure the air-cushioned face mask airtight^{12,14} while holding the airway open and the other hand to squeeze the self-inflating bag. The use of a BVM is a difficult technique to master and requires both hands of **the person performing the ventilation** as well as full attention to ensure the efficacy of ventilation. It is difficult to grasp the mask to keep the airway open with one hand while squeezing the bag with the other hand¹⁵. One hand is used to carefully squeeze the bag to provide adequate Vt^{13,16} and RR while being careful not to over pressurize thereby causing gastric insufflation or barotrauma. Because this task requires the entire attention of the person in charge of the ventilation, **a second person is needed to perform additional care tasks** such as administering medications.

Since it is difficult to provide a consistent RR during manual ventilation^{2,13}, bag-valve ventilation may have adverse clinical consequences^{7,10}. Alterations in the patient's respiratory status and arterial blood gas values can occur^{2,17}. The technique can result in **hyperventilation or hypoventilation**^{2,8,13,18,19}, and associated shifts in blood pH toward alkalosis or acidosis, respectively.

Dumont et al. findings indicate that hypocarbia due to hyperventilation¹⁸ and hypercarbia as a result of hypoventilation following **traumatic brain injury** (TBI)¹⁹⁻²¹ increase the risk of mortality; and normocarbia following TBI decreases the risk of in-hospital mortality. The authors propose that abnormal physiologic states of hypercarbia and hypocarbia may induce secondary injury by inducing suboptimal oxygenation of brain tissue during prehospital management. Because CO₂ arterial pressure might affect cerebral blood flow, both hypocapnia and hypercapnia should be avoided². Adjustment of ventilation to achieve

normocarbia can be facilitated by monitoring the end-tidal CO₂ during transport^{3,22}. Normocarbia appears to decrease the risk of mortality by maintaining a normal physiologic state and minimizing secondary brain injury. Therefore, normoventilation in prehospital setting is critical⁸.

In addition, the vast majority of patients with return of spontaneous circulation need assisted ventilation. Hypoxemia must be prevented, considering the evidence of harm after myocardial infarction and the possibility of increased neurological injury after **cardiac arrest**, it is recommended that the inspired oxygen concentration maintains arterial blood oxygen saturation in the range of 94–98%²². Moreover, during CPR the physiological consequences of hyperventilation⁸ and of prolonged ventilation intervals result in a persistently positive intrathoracic pressure during the decompression phase, thereby decreasing cardiac preload and cardiac output and impeding right ventricular function. Increased Vt is also known to adversely affect cardiac output⁹.

Accurate breath rates are difficult to maintain because delivering assisted ventilation, especially under frenzied environments⁸ and difficult clinical conditions, alters time perception¹³. Successful mask ventilation should give a Vt of 6-8 mL/Kg of ideal body weight⁴ and an upper limit pressure of 20-25 cmH₂O. Ensuring adequate ventilating pressures, to achieve therapeutic tidal volumes during manual ventilation, is also difficult. As mentioned before, **excessive pressures** can cause a decrease in cardiac preload^{3,9,10}, cardiac output and coronary perfusion in hemorrhagic shock and cardiac arrest; it can also cause secondary brain injury^{2,3,20,21} and gastric insufflation^{16,23,24}. Also, bag-valve ventilation may result in **lung injury**²⁵ suggesting an increased potential for **barotrauma**^{2,13}. However, insufficient pressures

can lead to inadequate ventilation and concomitant hypoxia and hypercapnia. Inappropriate mask ventilation in respiratory compromised patients is a medical emergency.

Tidal volumes must be large enough to assure oxygen delivery with positive-pressure ventilation; however, applying **high airway pressures** that exceed the lower esophageal sphincter pressure can force ventilating gas into the stomach¹⁶. **Gastric inflation** can elevate the diaphragm, restrict lung movement, and decrease respiratory system compliance⁴. It also increases the risk of regurgitation and aspiration of stomach contents¹³, particularly in patients who received concurrent chest compressions for CPR, victims of drowning, and cases in which lower esophageal sphincter pressure dropped severely due to prolonged hypoxia.

Mechanical ventilation before arrival to the emergency department

Due to the variability in respiratory parameters and the possible deterioration of blood gases and changes in hemodynamics during BVM use^{3,17}, mechanical ventilation with **portable devices** have become the preferred method for transporting patients in prehospital environments^{2,12,22}. In those situations, utilizing automated feedback control, even with noninvasive ventilation, is more efficacious and safer than manual ventilation¹³.

In 2020, Fogarty et al.¹³ explained that healthcare providers using emergency transport ventilators (ETV) control the inflation time, respiratory frequency and tidal volume limiting peak airway pressures, therefore can mitigate barotrauma and/or stomach inflation. **Prehospital noninvasive ventilation** (NIV) limits the pressures to safe ranges avoiding over pressurization and overdistension of the lungs. It utilizes positive end expiratory pressure (PEEP) during exhalation avoiding the collapse of diseased alveoli. Therefore, the main objective

of mechanical ventilation is to maximize gas exchange, oxygen delivery to the lungs and elimination of CO₂. The authors state that NIV improves alveolar ventilation while employing lung protective ventilation (< 20 cmH₂O).

In the prehospital setting the use of ETV should provide a stable minute ventilation and free up a caregiver to perform other tasks^{4,13}. By using these devices, EMS personnel can avoid exceeding the opening pressure of the lower esophageal sphincter and introducing air into the stomach. In addition, **ETV provide accurate and precise tidal volumes**. Also important, the use of these devices frees EMS providers to turn their attention to other patient care issues, and improve the overall patient care³. Prehospital use of ETV in patients with acute respiratory failure results in a **fast and significant improvement of vital functions**²⁶ such as blood pressure, heart rate, breathing rate, and oxygenation. This stabilizing effect seems to be largely independent of the duration of the ventilatory support.

O-Two Medical Technologies Inc. has overcome the challenges of manual ventilation by introducing the SMART Bag[®] MO Bag-Valve-Mask Resuscitator; however, this BVM controls flow assisting the rescuer to deliver proper RR and Vt. It does not eliminate the variability on the squeeze method given by the operator specially under hectic environments and difficult clinical conditions that affect ventilation parameters inadvertently.

Avoiding adverse events during transport is the primary goal in providing a safe transport. While many transports are uneventful, sometimes the patient's condition deteriorates as a consequence of the underlying disease progression³. Patient monitoring is difficult during transport, due to noise, limited space for the caregivers, and low light. ETV could resolve most of these issues¹⁷, portable ventilators let the rescuer focus on maintaining a mask seal,

observing the patient and checking live-monitoring parameters, without having to worry about squeezing a BVM correctly. **All patients requiring emergency ventilation must be adequately monitored**¹⁷, including continuous monitoring of end-tidal carbon dioxide concentrations²². During prolonged resuscitative efforts the use of ETV, pneumatically powered and time- or pressure-cycled, allows the EMS team to perform other tasks while providing adequate ventilation and oxygenation. Nonetheless, providers should always have a bag-mask device available for backup¹². Ideally, as mentioned before, SMART Bag[®] MO Bag-Valve-Mask offers better manual ventilation. It provides consistent controlled ventilations, flowrate and airway pressure while almost completely eliminating the risks associated with conventional BVM ventilation.

In Canada, McLeod and collaborators²⁷ reviewed the effects of prehospital NIV on in-hospital mortality. The authors conducted the study in adults with severe respiratory distress, due to acute cardiogenic pulmonary edema (CPE), exacerbation of chronic obstructive pulmonary disease (COPD), asthma, or pneumonia. They compared prehospital NIV and standard care and described that **prehospital noninvasive positive-pressure ventilation reduces in-hospital mortality** and the need for invasive ventilation. The authors recommend that continuous positive airway pressure (CPAP) should be considered as first line intervention, because it is cheap and easy to implement in clinical practice.

Recently, Mario Hensel and colleagues²⁶ in Germany reported in a study that **early initiation of NIV is important for successful treatment of acute respiratory failure**. Patients can benefit from initiation of non-invasive treatment already in the prehospital setting to avoid intubation and to improve patient outcomes. An important concern regarding the prehospital NIV is the

additional time required to establish the equipment and to adapt the NIV-mask to the face of the patient, when avoiding a significant leakage. However, the authors emphasized that the extra effort associated with a prolonged on-scene time is justified aiming to improve vital functions. Therefore, when the ambulance transport time is short, as in most metropolitan areas, **the prehospital application of NIV is justified even with short times from arrival at scene to hospital handover**.

Portable ventilators are becoming increasingly robust in capability, and their prehospital use is more widespread for out-of-hospital emergency transports. Portable mechanical ventilators are used in patients with acute respiratory failure or some respiratory compromise who need ventilatory support while transported. Around the world, prehospital use of ETV is largely variable, some countries limit their use to air medical transport and critical care transfer ambulances²⁸. Other countries mandate the use of prehospital ETV for patients who require ventilation providing advanced medical care.

In that regard, in 2001 Wayne et al.³ described concepts and application of prehospital ventilation and concluded: ***“providing a ventilator for every ambulance should be considered because alternate techniques for ventilation have proved inadequate.”***

As previously explained, **mechanical ventilators offer a safe alternative to manual ventilation** and allow healthcare providers to deliver consistent and controlled tidal volumes at a determined rate with or without additional medications. ETV should provide battery power sufficient to finish the transport; full ventilatory support; controlled Vt and RR; stable Vt with lung compliance changes; a disconnect alarm; airway pressure monitoring; stable PEEP; and FiO₂ up to 1.0. Other desirable features for transport ventilators are: lightweight, and easy to use³; low gas consumption; easy to trigger;

both volume and/or pressure modes; adjustable FiO₂; and ability to operate from a 50 psi O₂ source¹⁷.

Portable mechanical ventilators such as the **eSeries automatic transport ventilators** offered by O-Two Medical Technologies Inc.²⁹ provide controlled ventilation for resuscitation and patient transport. These electronically controlled, pneumatically powered devices give a range of ventilation solutions for all patient respiratory requirements during prehospital care. The eSeries ventilators are compact, lightweight, ideal for long transports with long battery operation time. They offer volume and/or pressure modes, continuous parameter monitoring, adjustable FiO₂ up to 1.0 and low gas consumption. eSeries transport ventilators ensure standardization of ventilation and adherence to international guidelines as they count with all necessary and desirable features for prehospital positive-pressure treatment.

In conclusion, ventilation assistance is an important aspect when managing out-of-hospital medical emergencies involving acute respiratory failure. Advantages of prehospital noninvasive positive-pressure ventilation are that it delivers consistent and controlled Vt and RR, it directly affects patients' potential for recovery, promote oxygenation, and may protect against aspiration; these benefits are difficult to achieve with BVM ventilation. In other words, optimal controlled ventilation is crucial since it allows consistent Vt, respiratory frequency, and limited airway pressures minimizing the risks of complications such as hyperventilation, barotrauma, gastric inflation and hypoventilation.

References

1. Jacobs P E, Grabinsky A. Advances in prehospital airway management. *Int J Crit Illn Inj Sci* 2014; 4(1): 57–64.

2. Auble T E, Menegazzi J J, Nicklas K A. Comparison of automated and manual ventilation in a prehospital pediatric model, *Prehosp Emerg Care* 1998; 2:2, 108-111.
3. Wayne M A, Delbridge T R, Ornato J P, et al. Concepts and application of prehospital ventilation. *Prehosp Emerg Care* 2001; 5:1, 73-78.
4. Cave D M, Gazmuri R F, Otto C W, et al. Part 7: CPR Techniques and Devices 2010 American Heart Association Guidelines for Cardiopulmonary Resuscitation and Emergency Cardiovascular Care. *Circulation* 2010; 122[suppl 3]:S720 –S728.
5. L'her E, Marjanovic N, Le floch S, et al. Automatic vs. manual bag valve resuscitation: An experimental bench test. (P2-85) Abstracts – 17th World Congress on Disaster and Emergency Medicine. *Prehosp Disaster Med*; 26:1, s164.
6. Lee H M, Cho K H, Choi Y H, et al. Can you deliver accurate tidal volume by manual resuscitator? *Emerg Med J* 2008; 25:632–634.
7. Nikolla D, Lewandowski T, Carlson J, et al. Mitigating hyperventilation during cardiopulmonary resuscitation. *Am J Emerg Med* 2016; 34:643-646.
8. Pitts S, Kellermann A L. Hyperventilation during cardiac arrest. *Lancet* 2004; 364(9431):313-5.
9. Aufderheide T P, Sigurdsson G, Pirralo R G, et al. Hyperventilation-induced hypotension during cardiopulmonary resuscitation. *Circulation* 2004; 109:1960-1965.
10. Aufderheide T P, Lurie K G. Death by hyperventilation: A common and life-threatening problem during cardiopulmonary resuscitation. *Crit Care Med* 2004; Suppl.32:9.
11. Neth M R, Benoit J L, Stolz U, et al. Ventilation in simulated out-of-hospital cardiac arrest resuscitation rarely meets guidelines. *Prehosp Emerg Care* 2020; 00:000–000.
12. Neumar R W, Otto C W, Link M S; et al. Part 8: Adult advanced cardiovascular life support 2010 American Heart Association guidelines for cardiopulmonary resuscitation

- and emergency cardiovascular care circulation 2010; 122[suppl 3]:S729–S767.
13. Fogarty M, Kuck K, Orr J, et al. A comparison of controlled ventilation with a noninvasive ventilator versus traditional mask ventilation. *J Clin Monit Comput* 2020; 34:771–777.
 14. Na J U, Han S K, Choi P C, et al. Influence of face mask design on bag–valve–mask ventilation performance: A randomized simulation study. *Acta Anaesthesiol Scand* 2013; 57:1186–1192.
 15. Kroll M, Das J, Siegler J. Can altering grip technique and bag size optimize volume delivered with bag-valve-mask by emergency medical service providers? *Prehosp Emerg Care* 2019; 23:2, 210-214.
 16. Fitz-Clarke, J R. Fast or slow rescue ventilations: A predictive model of gastric inflation. *Respir Care* 2018; 63(5):502–509.
 17. Blakeman T C, Branson R D. Inter- and Intra-hospital transport of the critically ill. *Respir Care* 2013; 58(6):1008 –1021.
 18. Davis D P, Dunford J V, Poste J C, et al. The Impact of hypoxia and hyperventilation on outcome after paramedic rapid sequence intubation of severely head-injured patients. *J Trauma* 2004; 57:1–10.
 19. Woods B, Ward S, Lindsell C J, et al. Mechanical ventilation of severe traumatic brain injury patients in the prehospital setting. *Air Med J* 2020; 39:410-413.
 20. Dumont T M, Visioni A J, Rughani A I, et al. Inappropriate prehospital ventilation in severe traumatic brain injury increases in-hospital mortality. *J Neurotrauma* 2010; 27:1233–1241.
 21. Warner K J, Cuschieri J, Copass M K, et al. The impact of prehospital ventilation on outcome after severe traumatic brain injury. *J Trauma* 2007; 62:1330 –1338.
 22. Hock Ong M E, Perkins G D, Cariou A. Out-of-hospital cardiac arrest 2. Out-of-hospital cardiac arrest: prehospital management. *Lancet* 2018; 391:980–88.
 23. Spoormans I, Van Hoorenbeeck K, Balliu L, et al. Gastric perforation after cardiopulmonary resuscitation: Review of the literature. *Resuscitation* 2010; 81:272-280.
 24. Smally A J, Ross M J, Huot C P. Gastric rupture following bag-valve-mask ventilation. *J Emerg Med* 2002; 22(1):27–29.
 25. Stephens R J, Siegler J E, Fuller B M. Mechanical ventilation in the prehospital and emergency department environment. *Respir Care* 2019; 64(5):595–603.
 26. Hensel M, Strunden M S, Tank S, et al. Prehospital non-invasive ventilation in acute respiratory failure is justified even if the distance to hospital is short. *Am J Emerg Med* 2019; 37:4, 651-656.
 27. Mal S, McLeod S, Iansavichene A, et al. Effect of out-of-hospital noninvasive positive-pressure support ventilation in adult patients with severe respiratory distress: A systematic review and meta-analysis. *Ann Emerg Med* 2014; 63:600-7.
 28. El Sayed M, Tamim H, Mailhac A, et al. Ventilator use by emergency medical services during 911 calls in the United States. *Am J Emerg Med* 2018; 36:5, 763-768.
 29. O_Two Automatic Electronic Ventilator. <https://otwo.com/products/automatic-transport-ventilators/automatic-transport-ventilators-2/eseries-ventilators/>. Accessed on September 2021.

O-Two Medical Technologies Inc.
 45A Armthorpe Road,
 Brampton, ON, Canada L6T 5M4
 Toll Free: (800) 387 3405
 Tel: (905) 792 OTWO (6896)
 Fax: (905) 799 1339
 Emails: resuscitation@otwo.com
 clinical@otwo.com
 Website: www.otwo.com