



O-Two Oxygen Flow Controllers

OXYGEN THERAPY TREATMENT FOR SINGLE PATIENT INCIDENTS

When you need to provide oxygen therapy the O-Two adjustable oxygen therapy flow controller provides high-quality manufacturing and performance.

The O-Two adjustable oxygen therapy flow controller is a calibrated, 12-position, fixed orifice, non-gravity sensitive flow controller. The flow controller is available in a complete range of Chemetron, Ohio, female swivel DISS or 1/8 NPT female pipe inlet connections and DISS fitting valve outlet connections.

The flow controller is a lightweight, unibody, ergonomically designed assembly that can accommodate humidifiers. Ideal for all applications where adjustable oxygen therapy is to be provided, the O-Two adjustable oxygen therapy flow controller is unbeatable in quality, size, and performance

- OUTSTANDING IMPACT RESISTANT
- ALL-METAL CONSTRUCTION
- LARGE VIEWING WINDOW
- RIBBED CONTROL KNOB
- POSITIVE 12 CLICK-STOP POSITIONS



SPECIFICATIONS

THERAPY FLOW OPTION 1	0, 0.5, 1, 2, 3, 4, 5, 6, 8, 10, 12, 15 (+/- 10%)
THERAPY FLOW OPTION 2	0, 1, 2, 3, 4, 6, 8, 10, 12, 15, 20, 25 (+/- 10%)
MAXIMUM RATED INPUT PRESSURE	50 PSI
OPERATING TEMPERATURE RANGE	-18°C to 50°C (0°F to 122°F)
STORAGE TEMPERATURE RANGE	-40°C to 60°C (-40°F to 140°F)
SIZE (EXCLUDING OUTLET AND INLET FITTING) (H X W)	2.0 x 1.37 in / 5.1 x 3.5 mm
WEIGHT	4.0 oz / 0.1 kg

ORDERING INFORMATION

01FC0115D	Adjustable 0 - 15 L/Min Click Stop Flow Controller	Each
01FC0125D	Adjustable 0 - 25 L/Min Click Stop Flow Controller	Each

THIS PRODUCT HAS A TWO YEAR WARRANTY AGAINST MANUFACTURERS DEFECTS.



www.otwo.com

45A Armthorpe Road, Brampton, ON, Canada, L6T 5M4
 Telephone: +1 905 792-OTWO (6896) N.A. Toll Free: +1 800 387 3405
 Facsimile: +1 905 799 1339 Email: resuscitation@otwo.com