




Rescue Breather CPR Protective Barrier

Bacterial and viral protection
during CPR Training

The single-use O-Two Rescue Breather CPR protective barrier is ready when needed. Designed as a first-response device for use during CPR, the O-Two Rescue Breather provides a physical barrier between the rescuer and the patient. This single-use, cost-effective quality device provides rescuer protection by using a one-way flapper valve and filter to reduce the risk of cross-contamination.

Compact and lightweight, the O-Two Rescue Breather can be easily stored in first aid kits, purses, sports bags, glove compartments, and keychain pouches. Packaged complete with instructions in a clear polyethylene Ziploc® bag.



One size fits
All patients

Ergonomic patient
mouthpiece

Round rescuer
mouthpiece

One-way valve
with filter

Large faceshield gives excellent coverage yet
allows the rescuer to see the patient's face

ORDERING INFORMATION

01RB1000

O-Two rescue breather in ziplock plastic bag

Case/50

NOTE: This device is for emergency use only by persons trained in CPR and in the use of this device. This device does not guarantee 100% protection against disease transmission

www.otwo.com

45A Armthorpe Road, Brampton, ON, Canada, L6T 5M4

Telephone: +1 905 792-OTWO (6896) N.A. Toll Free: +1 800 387 3405

Facsimile: +1 905 799 1339 Email: resuscitation@otwo.com