

Guedel Airway Kit

The O-Two Guedel Airway Kit comes in a practical acrylic box containing 8 ergonomic units. Oropharyngeal airways come in various sizes, from infant to adult, and are commonly used in pre-hospital emergency care and for short-term airway management post-anesthetic or when manual methods are inadequate to maintain an open airway.



Choose the proper size by measuring from the first incisors to the angle of the jaw.



The oral airway is then inserted into the person's mouth upside down.



Once contact is made with the back of the throat, the airway is rotated 180 degrees,



allowing for easy insertion, and assuring that the tongue is secured.

8 Different sizes
for different ages

Safe and
ergonomic

Practical
acrylic box

ORDERING INFORMATION

01AM3000

O-Two Guedel Airway Kit

Each

www.otwo.com

45A Armthorpe Road, Brampton, ON, Canada, L6T 5M4

Telephone: +1 905 792-OTWO (6896) N.A. Toll Free: +1 800 387 3405

Facsimile: +1 905 799 1339 Email: resuscitation@otwo.com