



O-Two CPR Personal Protection Devices

Bacterial and viral protection
during CPR Training

O-Two CPR Bio-Barrier® Training Faceshields

The two-way biological filter of the O-Two CPR Bio-Barrier training faceshield will prevent the passage of over 99% of the viruses and bacteria that cause diseases such as tuberculosis, hepatitis and HIV*. It also acts as a physical barrier to moisture, helping to prevent the transfer of body fluids.

During CPR training, the large, pliable — yet extremely strong — faceshield with graphic illustrations assists the student in correctly positioning the O-Two CPR Bio-Barrier training faceshield.

O-Two CPR Bio-Barrier training faceshields are packaged on a roll, allowing each student to tear off her or his own personal protective device.

*Independent laboratory tested



Features

Bacterial/viral filtration

Clear, oversized faceshield

Faceshield is pliable yet strong

Positioning illustrations printed on the faceshield

Ordering Information

01RB5050 CPR Bio-Barrier Training Faceshield

Dispenser roll of 50 (4 rolls per case)



45A Armthorpe Road, Brampton, ON, Canada, L6T 5M4

Telephone: +1 905 792-OTWO (6896)

N.A. Toll Free: +1 800 387 3405 Facsimile: +1 905 799 1339

Web: otwo.com Email: resuscitation@otwo.com



CPR Bio-Barrier is a registered trademark of O-Two Medical Technologies Inc.

15PL7000-A Rev. 110_16