



O-Two Multi Patient Oxygen Therapy Unit and Flow Controls

Oxygen Therapy Treatment for Single Patient and
Multiple Casualty Incidents

O-Two Multi Patient Unit

Oxygen Therapy Treatment for Multiple Casualty Incidents

When multiple casualty incidents occur, EMS providers often find their resources stretched in terms of equipment as well as personnel. For incidents where a number of patients require oxygen therapy simultaneously, the O-Two multi patient unit meets the challenge.

With six 0-15 L/min adjustable therapy flow controls and the capacity to operate an automatic resuscitator or demand valve resuscitator from the 9/16 in DISS self-sealing outlet on the end of the unit, the flexibility offered by the compact O-Two multi patient unit is second to none.

The O-Two multi patient unit attaches to any oxygen regulator with a 9/16 in DISS outlet that supplies 50 PSI pressure, with a minimum output flow rate of at least 160 L/min.



Features

All metal construction

Six individual adjustable therapy flow controls each providing a 0-15 L/min flow rate

One 9/16 in DISS self-sealing outlet for the attachment of an automatic ventilator or demand valve

Twenty-five foot supply hose is standard, allowing the unit to be placed where it is needed

Supplied with six adult medium concentration oxygen therapy facemasks with sure-flow tubing

Specifications

Input Pressure:	50 PSI
Operating Temperature:	-18 to 50°C (0 to 122°F)
Storage Temperature:	-40 to 60°C (-40 to 140°F)
Relative Humidity for Storage and Operating Use:	15 to 95%
Input Connection:	9/16 in DISS
Output Connection:	9/16 in DISS self-sealing outlet
Therapy Flow Connections:	Firtree
Therapy Flow Range:	0, 0.5, 1, 2, 3, 4, 5, 6, 8, 10, 12, 15 L/min
Weight (Multi Patient Unit only):	19.2 oz / 0.45 kg (approx.)
Dimensions (inches):	8.0 x 5.0 x 2.0 in (approx.)
(millimeters):	203 x 127 x 50 mm (approx.)

Ordering Information

01OR6304 Multi Patient Unit

Incl., Twenty-five foot supply hose, six (6) 0-15 L/min therapy flow controls, one (1) self-sealing 9/16 in DISS outlet, six (6) medium concentration masks with tubing supplied in a soft-pack carrying case.



45A Armthorpe Road, Brampton, ON, Canada, L6T 5M4

Telephone: +1 905 792-OTWO (6896)

N.A. Toll Free: +1 800 387 3405 Facsimile: +1 905 799 1339

Web: otwo.com Email: resuscitation@otwo.com



15PL1074 Rev. 10 10_16