



O-Two CPR Personal Protection Devices

Bacterial and viral protection during training and CPR

O-Two Rescue Breather® CPR Protective Barrier

The single-use O-Two Rescue Breather CPR protective barrier is ready when you need it. Designed as a first-response device for use during CPR, the O-Two Rescue Breather provides a physical barrier between the rescuer and the patient. This single-use, cost-effective quality device provides rescuer protection by the use of a one-way flapper valve and filter to reduce the risk of cross-contamination.

Compact and lightweight, the O-Two Rescue Breather can be easily stored in first aid kits, purses, sports bags and glove compartments. Packaged complete with instructions in a clear polyethylene Ziploc® bag, the O-Two Rescue Breather is also available in a keychain pouch or a key chain/belt pouch complete with a pair of gloves, ensuring that you will always carry the level of protection you deserve when helping to save a life.



Features

- One size fits all patients
- Ergonomic patient mouthpiece
- Round rescuer mouthpiece
- One-way valve with filter
- Large faceshield gives excellent coverage yet allows the rescuer to see the patient's face

Ordering Information

- 01RB1000-cs Rescue Breather CPR Device**
Clear polyethylene bag (Case of 50)
- 01RB2010 Red Rescue Breather CPR Device**
Red nylon keychain pouch (Case of 50)
- 01RB2010 Black Rescue Breather CPR Device**
Black nylon keychain pouch (Case of 50)
- 01RB2020-cs Rescue Breather CPR Device**
Nylon keychain/belt pouch with gloves (Case of 50)
(Available in Red, Black and Blue)



45A Armthorpe road, Brampton, ON, Canada, L6T 5M4
Telephone +1 905 792-OTWO (6896) N.A. Toll Free +1 800 387 3405
Facsimile +1 905 799 1339 Web otwo.com
Email resuscitation@otwo.com

Rescue Breather is a registered trademark of O-Two Medical Technologies Inc.
Ziploc is a registered trademark of SC Johnson & Son, Inc.



15PL7000-1 Rev. 8/11/16