



# CPR Bio-Barrier® Training Faceshields

# BACTERIAL AND VIRAL PROTECTION DURING CPR TRAINING

The two-way biological filter of the O-Two CPR BioBarrier training faceshield will prevent the passage of over 99% of the viruses and bacteria that cause diseases such as tuberculosis, hepatitis and HIV\*. It also acts as a physical barrier to moisture, helping to prevent the transfer of body fluids.

During CPR training, the large, pliable — yet extremely strong — faceshield with graphic illustrations assists the student in correctly positioning the O-Two CPR BioBarrier training faceshield.

O-Two CPR Bio-Barrier training faceshields are packaged on a roll, allowing each student to tear off her or his own personal protective device.



BACTERIAL/VIRAL FILTRATION

CLEAR, OVERSIZED FACESHIELD

FACESHIELD IS PLIABLE YET STRONG

POSITIONING ILLUSTRATIONS PRINTED ON THE FACESHIELD

## ORDERING INFORMATION

**01RB5050** - Bio-Barrier Training Faceshield  
Dispenser roll of 50 (4 rolls per case)

**o\_two** controlled ventilation™

[www.otwo.com](http://www.otwo.com)

rpe Road, Brampton, ON, Canada, L6T 5M4

Telephone: +1 905 792-OTWO (6896) N.A. Toll Free: +1 800 387 3405

Facsimile: +1 905 799 1339 Email: [resuscitation@otwo.com](mailto:resuscitation@otwo.com)

