

# Rescue Breather CPR Protective Barrier



# BACTERIAL AND VIRAL PROTECTION DURING CPR TRAINING

The single-use O-Two Rescue Breather CPR protective barrier is ready when you need it. Designed as a first-response device for use during CPR, the O-Two Rescue Breather provides a physical barrier between the rescuer and the patient. This single-use, cost-effective quality device provides rescuer protection by the use of a one-way flapper valve and filter to reduce the risk of crosscontamination.

Compact and lightweight, the O-Two Rescue Breather can be easily stored in first aid kits, purses, sports bags and glove compartments. Packaged complete with instructions in a clear polyethylene Ziploc® bag, the O-Two Rescue Breather is also available in a keychain pouch or a key chain/belt pouch complete with a pair of gloves, ensuring that you will always carry the level of protection you deserve when helping to save a life.



ONE SIZE FITS ALL PATIENTS

ERGONOMIC PATIENT MOUTHPIECE

ROUND RESCUER MOUTHPIECE

ONE-WAY VALVE WITH FILTER

LARGE FACESHIELD GIVES EXCELLENT COVERAGE YET  
ALLOWS THE RESCUER TO SEE THE PATIENT'S FACE

## ORDERING INFORMATION

**01RB1000-CS** - Rescue Breather CPR Device

[Clear polyethylene bag \(Case of 50\)](#)

**01RB2010** - Red Rescue Breather CPR Device

[Red nylon keychain pouch \(Case of 50\)](#)

**01RB2010** - Black Rescue Breather CPR Device

[Black nylon keychain pouch \(Case of 50\)](#)

**01RB2020-CS** - Rescue Breather CPR Device

[Nylon keychain/belt pouch with gloves \(Case of 50\)](#)

(Available in Red, Black and Blue)

**o\_two** controlled  
ventilation™

[www.otwo.com](http://www.otwo.com)

rpe Road, Brampton, ON, Canada, L6T 5M4

Telephone: +1 905 792-OTWO (6896) N.A. Toll Free: +1 800 387 3405

Facsimile: +1 905 799 1339 Email: [resuscitation@otwo.com](mailto:resuscitation@otwo.com)

CE 0120

15PL7000-1 Rev. 9 06\_18