

Resuscitation and Anesthesia Facemasks



EFFECTIVE AND ECONOMICAL

With a wide choice in style, material and size, O-Two provides a comprehensive range of resuscitation and anesthesia facemasks to meet any respiratory requirement.

From the truly “one-size-fits-all” Universal Facemask (available in both PVC and silicone) to the disposable, inflated-cuff range for neonates to adults, the high quality and superior design of these masks makes them the perfect choice for both hospital and pre-hospital use.

The silicone range includes donut-style neonate/infant masks as well as the flanged design of the O-Two Everseal Edge range. These are ideal anesthesia masks as they withstand repeated autoclaving.

The three sizes of clear PVC Everseal masks are ideal for pre-hospital use where disinfection after use is normally undertaken by cold chemical sterilant such as an activated glutaraldehyde or a 10% bleach solution.

Where a single-use mask is required, the disposable, inflated cuff range, meets the user’s needs. These masks come with clear plastic bodies and soft sealing cushions, which reduce the chance of cuff failure. The inflation pressure of the cuff is calculated to provide a superior seal when hand pressure is applied to the mask. Our single-use cuffed masks are provided with a head harness hook ring.

All masks have standard 22 mm or 15 mm (infant/neonate) inlet connectors and come conveniently packaged in cases of 12.



ORDERING INFORMATION

02FM4999-cs - Universal Mask

[PVC \(Case of 12\)](#)

02FM5000-cs - Medium Mask (#2)

[With anatomical chin cup \(PVC\) \(Case of 12\)](#)

02FM5001-cs - Large Mask (#3)

[With anatomical chin cup \(PVC\) \(Case of 12\)](#)

02FM5097-cs - Universal Mask

[Silicone \(Case of 12\)](#)

02FM5098-cs - Medium Mask (#2)

[with anatomical chin cup \(Silicone\) \(Case of 12\)](#)

02FM5099-cs - Large Mask (#3)

[with anatomical chin cup \(Silicone\) \(Case of 12\)](#)

02FM5100-cs - Toddler/Child Mask

[Silicone \(Case of 12\)](#)

02FM5101-cs - Infant Mask

[Silicone \(Case of 12\)](#)

02FM5102-cs - Neonate Mask

[Silicone \(Case of 12\)](#)

02FM5203-Cs - Mask Disposable Cuffed Child #2

[with Hook Ring \(Case of 12\)](#)

02FM5205-Cs - Mask Disposable Cuffed Infant #1

[with Hook Ring \(Case of 12\)](#)

02FM5206-Cs - Mask Disposable Cuffed Neonate #0

[with Hook Ring \(Case of 12\)](#)

02FM5207-Cs - Mask Disposable Cuffed Child #3

[with Hook Ring \(Case of 12\)](#)

02FM5208-Cs - Mask Disposable Cuffed Lge Adult #5

[with Hook Ring \(Case of 12\)](#)

02FM5210-Cs - Mask Disposable Cuffed Adult #4

[with Hook Ring \(Case of 12\)](#)

o_two controlled
ventilation™

www.otwo.com

45A Armthorpe Road, Brampton, ON, Canada, L6T 5M4

Telephone: +1 905 792-OTWO (6896) N.A. Toll Free: +1 800 387 3405

Facsimile: +1 905 799 1339 Email: resuscitation@otwo.com

