



# O-Two Oxygen Therapy Kits

**o\_two** controlled  
ventilation™

# OXYGEN THERAPY ON THE MOVE!

O-Two Medical Technologies Oxygen Therapy Kits are ready for any breathing emergency when seconds count. They can be ordered either in a Emergency Responder Moulded Plastic Carrying Case constructed in high visibility International Orange or our Royal Blue Softpack Clamshell Carrying Case with Sling Style Shoulder Strap for "hands-free" operation.



The O-Two Medical Technologies Adjustable Oxygen Therapy Regulators feature high quality manufacturing and performance. The "Aluminum body / Brass Core" construction improves product safety by reducing the risk of exacerbating regulator fires should they occur.



The Emergency Responder Moulded Plastic Carrying Case has a built in carrying handle that allows the kit to be carried in a horizontal position and remains in a stable position during use. A moulded transport hole in the end of the case allows the oxygen therapy mask to be used with the case lid closed during patient transport. The case design makes the kit lightweight and buoyant should the kit fall into the water during an emergency. It will stay afloat for a short period of time allowing for immediate retrieval.



The Softpack Clamshell Carrying Case is constructed with extra thick closed cell foam, elastic airway holders, double reinforced abrasion resistant nylon bottom protecting your equipment from the rigors of Emergency use. It features a large two way, self-repairing zipper that allow the case to be opened fully during any breathing emergency. Large, oversized, interior pockets allow the easy access storage of oxygen therapy disposables and accessories. The outside features reflective material for visibility and double carrying handles for security and an adjustable shoulder strap sling that frees the rescuer's hands for other tasks.

THESE PRODUCTS ARE IDEAL FOR ALL APPLICATIONS WHERE ADJUSTABLE OXYGEN THERAPY FROM 0 - 25 L/MIN FLOWRATES ARE REQUIRED.

## SPECIFICATIONS

### MOULDED CARRY CASE

SIZE (APPROX):	L 23" x H 6.25" x W 10.25" (approx)
COLOR:	International Orange

### CLAMSHELL SOFTPACK

SIZE (APPROX):	L 19" x Dia. 8" (approx)
COLOR:	Royal Blue

### REGULATOR

MAX NORMAL INLET PRESSURE (P1):	3000 PSI (206.8 BAR)
OPERATING INLET PRESSURE:	500-3000 PSI (34.5-206.8 BAR)
NOMINAL OUTLET PRESSURE:	60 PSI (4.14 BAR)
MINIMUM DYNAMIC PRESSURE:	46 PSI
GAUGE ACCURACY:	2% of full scale + 8% of reading
STORAGE TEMPERATURE:	-28 to 60°C (-18.4 to 140°F)
OPERATING TEMPERATURE:	-18 to 50°C (-0.4 to 122°F)
SIZE (INCL., GAUGE):	6"L x 3"H x 1.5"W (approx.)
INCHES:	6"L x 3"H x 1.5"W (approx.)
MILIMETER:	152 x 76 x 38 (approx.)
WEIGHT (APPROX.):	27 oz / 748 g

## ORDERING INFORMATION

**01OT3000** - O-Two Medical Inhalator Kit c/w Adjustable 0 - 25 L/min Therapy Flow Regulator, Cylinder Wrench, Medium Concentration Facemask with Tubing in a Moulded Carrying Case (Cylinder not Included)

**01OT3200** - O-Two Medical Inhalator Kit c/w Adjustable 0 - 25 L/min Therapy Flow Regulator, Cylinder Wrench, Medium Concentration Facemask with Tubing in a Softpack Clamshell Carrying Case (Cylinder not included)

**01OT3300** - O-Two Medical Technologies Emergency Responder Kit c/w Adjustable 0-25 L/min Regulator, Disposable Easy Grip® Adult Resuscitator, Disposable High Concentration Adult and Child Oxygen Masks with Tubing, 2 Pairs Gloves in Softpack Clamshell Carrying Case with "D" Aluminum Oxygen Cylinder



[www.otwo.com](http://www.otwo.com)

rpe Road, Brampton, ON, Canada, L6T 5M4

Telephone: +1 905 792-OTWO (6896) N.A. Toll Free: +1 800 387 3405

Facsimile: +1 905 799 1339 Email: [resuscitation@otwo.com](mailto:resuscitation@otwo.com)

