



# Rescue Breather CPR Pocket Ventilator



## BACTERIAL AND VIRAL PROTECTION DURING TRAINING AND CPR

The O-Two Rescue Breather CPR pocket ventilator provides the rescuer with a safe and effective means to perform life-saving expired air resuscitation while offering maximum protection from the potential risks of disease transfer for both the patient and the rescuer.

The O-Two Rescue Breather CPR pocket ventilator comes standard with a low-resistance disposable one-way valve, with bacterial/viral filter, oxygen inlet, head-strap and carrying case. The filter ensures optimum bacterial protection, yet minimum airway resistance. The seamless air cushion mask eliminates seam failures and air loss during storage.

The compact size of the O-Two Rescue Breather CPR pocket ventilator allows it to fit into your pocket, purse, brief case, sports bag or glove compartment of your car — ready for any breathing emergency.

- TRANSPARENT UNIVERSAL CUSHION
- MASK FOR ADULTS TO INFANTS
- REPLACEABLE, CLEAR, ONE-WAY VALVE WITH BACTERIAL/VIRAL FILTER
- OXYGEN INLET WITH PROTECTIVE CAP
- HEAD-STRAP
- COMPACT, CLAMSHELL CASE WITH HOLD TABS

## ORDERING INFORMATION

**01PV1500W-CS** - Rescue Breather® CPR Pocket Ventilator  
Includes head-strap, oxygen inlet, one-way valve and bacterial/viral filter in a clamshell hard case (Case of 12)

**01PV1505W-CS** - Replacement One-Way Valves  
(Package of 10) (Also available in an unprinted clamshell case - **01PV1503W-CS**)



[www.otwo.com](http://www.otwo.com)

45A Armthorpe Road, Brampton, ON, Canada, L6T 5M4

Telephone: +1 905 792-OTWO (6896) N.A. Toll Free: +1 800 387 3405

Facsimile: +1 905 799 1339 Email: [resuscitation@otwo.com](mailto:resuscitation@otwo.com)

