



O-Two Oxygen Therapy Kits

“Oxygen therapy On The Move”!

O-Two Medical Technologies Oxygen Therapy Kits are ready for any breathing emergency when seconds count. They can be ordered either in a Emergency Responder Moulded Plastic Carrying Case constructed in high visibility International Orange or our Royal Blue Softpack Clamshell Carrying Case with Sling Style Shoulder Strap for “hands-free” operation.

The O-Two Medical Technologies Adjustable Oxygen Therapy Regulators feature high quality manufacturing and performance. The “Aluminum body / Brass Core” construction improves product safety by reducing the risk of exacerbating regulator fires should they occur.

These products are ideal for all applications where adjustable oxygen therapy from 0 - 25 L/min flowrates are required.



SPECIFICATIONS

Carrying Cases

Moulded Carry Case:

Size (approx): L 23" x H 6.25" x W 10.25" (approx)
Color: International Orange

Clamshell Softpack:

Size (approx): L 19" x Dia. 8" (approx)
Color: Royal Blue

Regulator:

Max Normal Inlet

Pressure (p1): 3000 PSI (206.8 BAR)

Operating Inlet

Pressure: 500-3000 PSI (34.5-206.8 BAR)

Nominal Outlet

Pressure: 60 PSI (4.14 BAR)

Minimum Dynamic

Pressure: 46 PSI

Gauge Accuracy: 2% of full scale + 8% of reading

Storage Temperature: -28 to 60°C (-18.4 to 140°F)

Operating

Temperature: -18 to 50°C (-0.4 to 122°F)

Size (incl., gauge): L x H x W

(inches): 6 x 3 x 1.5 (approx.)

(millimeters): 152 x 76 x 38 (approx.)

Weight (approx.): 27 oz / 748 g

The **Emergency Responder Moulded Plastic Carrying Case** has a built in carrying handle that allows the kit to be carried in a horizontal position and remains in a stable position during use. A moulded transport hole in the end of the case allows the oxygen therapy mask to be used with the case lid closed during patient transport. The case design makes the kit lightweight and buoyant should the kit fall into the water during an emergency. It will stay afloat for a short period of time allowing for immediate retrieval.

The **Softpack Clamshell Carrying Case** is constructed with extra thick closed cell foam, elastic airway holders, double reinforced abrasion resistant nylon bottom protecting your equipment from the rigors of Emergency use. It features a large two way, self-repairing zipper that allow the case to be opened fully during any breathing emergency. Large, oversized, interior pockets allow the easy access storage of oxygen therapy disposables and accessories. The outside features reflective material for visibility and double carrying handles for security and an adjustable shoulder strap sling that frees the rescuer's hands for other tasks.

ORDERING INFORMATION:

- 010T3000 O-Two Medical Inhalator Kit**
c/w Adjustable 0 - 25 L/min Therapy Flow Regulator, Cylinder Wrench, Medium Concentration Facemask with Tubing in a Moulded Carrying Case (Cylinder not Included)
- 010T3200 O-Two Medical Inhalator Kit**
c/w Adjustable 0 - 25 L/min Therapy Flow Regulator, Cylinder Wrench, Medium Concentration Facemask with Tubing in a Softpack Clamshell Carrying Case (Cylinder not included)
- 010T3300 O-Two Medical Technologies Emergency Responder Kit**
c/w Adjustable 0-25 L/min Regulator, Disposable Easy Grip® Adult Resuscitator, Disposable High Concentration Adult and Child Oxygen Masks with Tubing, 2 Pairs Gloves in Softpack Clamshell Carrying Case with “D” Aluminum Oxygen Cylinder



45A Armthorpe Road, Brampton, ON, Canada, L6T 5M4

Telephone: +1 905 792-OTWO (6896)

N.A. Toll Free: +1 800 387 3405 Facsimile: +1 905 799 1339

Web: otwo.com Email: resuscitation@otwo.com



Part #: 15PL7020 Rev. 4 10_16