



2024

EUROPEAN
CATALOG

o_two controlled[™]
ventilation

Effective January 1st 2024

CONTENTS

AUTOMATIC VENTILATION	04
ANALGESIA	09
IMMEDIATE CARE	12
OXYGEN DEMAND VALVES & KITS	18
PRESSURE REGULATORS & OXYGEN ADMINISTRATION	21
MOUNTING AND CARRYING SOLUTION	24

O-TWO MEDICAL TECHNOLOGIES INC. CUSTOMER SERVICE POLICY

1. CONDITIONS OF SALE:

- All orders are subject to acceptance by O-Two Medical Technologies Inc.
- New accounts require credit approval prior to shipment of goods.
- Minimum order value \$250.00
- Orders below this amount must be prepaid by wire transfer, bank draft Prices are subject to change without notice.
- O-Two reserves the right to cancel any order at any time, in its sole and absolute discretion, without liability for additional compensation to the buyer.
- The submission of a Purchase Order or payment of a Proforma Invoice constitutes the purchaser's acceptance of O-Two's customer service policy.

2. FREIGHT POLICY:

- Orders will be shipped collect via carrier of your choice or shipping costs may be prepaid by O-Two and applied to your invoice.
- All prices are Ex Works Brampton, Ontario.
- The shipment insurance is the purchaser's responsibility in all shipping terms.

3. PLACING ORDERS:

- Orders must be placed by e-mail or fax direct with O-Two Customer Service.
- All sales are in full case quantities only.
- Standard shipping is 2 weeks (unless otherwise notified by Customer Service).

4. BACKORDER POLICY:

- All backordered items will be shipped immediately upon availability, unless otherwise specified.

5. PAYMENT TERMS:

- Net 30 from date of invoice on approved credit.
- All other payments on non-credit accounts may be made by wire transfer bank draft prior to acceptance of any order.
- A 2% monthly interest charge on past due balances may be applied.

6. RETURN GOODS POLICY:

- A Return Authorization Number (RAN) must be obtained prior to returning any merchandise. Any product returned without a RAN will be refused. Merchandise shipped from O-Two Medical Technologies Inc. more than 30 days prior to return authorization request may not be returned.
- Freight charges on returned goods must be prepaid.
- Returned goods are subject to our inspection before credit is issued. If saleable merchandise is returned for credit there will be a 20% restocking charge deducted from the credit note.
- Where merchandise is returned as a result of O-Two Medical Technologies Inc. error or manufacturing defect, full credit or replacement of merchandise, including freight, at the discretion of O-Two Medical Technologies Inc. will be issued.
- The Authorization form with the Returned Goods Authorization number must be completed and enclosed with all returned goods along with a packing slip.
- Obsolete products or products made for special order cannot be canceled once order confirmation is sent by customer service, and are not returnable.

Bio Barrier®, Rescue Breather®, Equinox®, Statvac®, SMART BAG®, Easy Grip®, Revive-Aid® and CAREvent® are all registered trademarks of O-Two Medical Technologies Inc.

AUTOMATIC VENTILATION

eSeries Ventilators



01EVE500

e500 - Electronic Automatic Transport Ventilator c/w Disposable Patient Circuit, Resuscitation Mask, 6 Foot Supply Hose with 9/16" DISS Fittings, Calibrated Test Lung and Power Supply.

Each



01EVE700

e700 - Electronic Automatic Transport Ventilator c/w Disposable Patient Circuit, Resuscitation Mask, 6 Foot Supply Hose with 9/16" DISS Fittings, Calibrated Test Lung and Power Supply.

Each

eSeries Ventilators Accessories



PART N°	DESCRIPTION	UNITS
01CV8030-CS	O-Two Medical Single-Use Electronic Transport Ventilator 6 Foot Circuit with Protective Sleeve. For e500, e600 and e700	Case/10



01CV7035	O-Two Medical “e” Series Ventilator “Smart Mount” Multi-configuration Mounting Bracket (Ambulance Cot, Hospital Stretcher, Bed, Roll Stand)	Each
-----------------	---	------



01CV8040-CS	O-Two Medical “e” Series Ventilator Replacement Intake Filter/Cover	Case/10
--------------------	---	---------



01TA1852	O-Two Medical “e” Series Ventilator Replacement 1 Litre Test Lung with Compliance Restrictor	Each
-----------------	--	------



01FV4303-DISS	O-Two 6 Foot (1.85 Meter) O2 Supply Hose with 9/16 DISS Nut and 9/16” DISS Nut Ventilator Connection	Each
----------------------	--	------



01CV9100	O-Two Medical “e” Series Replacement Lithium Ion Replacement Battery	Each
-----------------	--	------



01CV0105	Power Supply - eSeries	Each
-----------------	------------------------	------



01CV0106	Power supply cord for eSeries power supply	Each
-----------------	--	------



01CV0102-EU	Power supply cord for eSeries power supply	Each
--------------------	--	------



01CV7050	eSeries Automatic Transport Ventilator Carrying Case - With sling-style shoulder strap (specifically designed for eSeries)	Each
-----------------	--	------

*Refer to page 34 to confirm the hose required for your country

CAREvent® Ventilators



PART N°	DESCRIPTION	UNITS
01CV6000	CAREvent® ATV+ Automatic Transport Ventilator c/w Adjustable 20 - 60 cm H ₂ O Pressure Relief System, Low Pressure Visual and Audible Alarms, Air Mix, Adjustable CPAP/PEEP, Breathing System Integrity (BSI) Alarm, 6' Supply Hose*, Deluxe Disposable Ventilation Circuit, Calibrated Test Lung and Universal Resuscitation Mask	Each
01CV7000	CAREvent® MRI Conditional Transport Ventilator c/w Adjustable 20 - 60 cm H ₂ O Pressure Relief System, Low Pressure Visual and Audible Alarm, Air Mix, Adjustable CPAP/PEEP, Breathing System Integrity (BSI) Alarm, 6' Supply Hose*, Deluxe Disposable Ventilation Circuit, Calibrated Test Lung & Universal Mask	Each
01CV3000	CAREvent® ALS Automatic Transport Resuscitator complete with Gurney Clip, 6 Foot Supply Hose* and Universal Resuscitation Mask	Each
01CV3000-CA	CAREvent® CA (Chemical Agent) Handheld Automatic Rescue Resuscitator complete with 6 Foot Supply Hose* and Universal Resuscitation Face Mask	Each

CAREvent® Accessories



PART N°	DESCRIPTION	UNITS
01CV8015-CS	O-Two Deluxe Disposable Ventilation Circuit with PEEP Port & anti-entrainment flapper (for use with ALS, ALS+, ATV+, MRI Compatible Ventilators)	Case/10
01CV8016-CS	O-Two CAREvent® MRI Ventilation Circuit with PEEP Port and 10 Foot Hose	Case/10
01CV8010-CS	O-Two Deluxe Disposable Ventilation Circuit with PEEP Port (without anti-entrainment flapper)	Case/10
10-55330-CS	Disposable PEEP Valve 5 - 20 cmH ₂ O	Case/12
17MP8690-CS	Air Intake Filter and Cover for CAREvent® ATV, ATV+, MRI Compatible Ventilators	Case/10
01PV1121	Duckbill diaphragm	Each
17MP1535-CS	Anti-entrainment ring	Case/10
17MP9039	Non-rebreathing valve cap assembly	Each
01FG6502	CAREvent® Handheld Automatic Time/Volume Cycled Resuscitator Head Harness System	Each












*Refer to page 34 to confirm the hose required for your country

ANALGESIA




PART N°	DESCRIPTION	UNITS
01EQ3000	eAdvantage® Adjustable Analgesic Gas Mixing and Delivery System with Supply Hoses (c/w DISS O ₂ and N ₂ O), Single-Use Patient Circuit with Scavenger Hose and Adapter with Mouthpiece	Each
01EQ1000	Equinox® Relieve N ₂ O/O ₂ Analgesic Gas Mixing and Delivery System c/w Supply Hoses (c/w DISS O ₂ and N ₂ O), Single-Use Patient Circuit with Scavenger Hose and Adapter with Mouthpiece	Each
01EQ7000E	Equinox® II Administration System c/w Regulator, Bleed Pin, 7 Foot Supply Hose, Demand Valve, Universal Face Mask, Cylinder Wrench	Case/10
01EQ7001E	Equinox® II Demand Valve c/w 7 Foot Supply Hose, Universal Resuscitation Face Mask	Case/10

Accessories

PART N°	DESCRIPTION	UNITS
 01EQ7006	Equinox® II Demand Valve Bleed Pin with Chain	Each
 17MP1528	Equinox® II Diaphragm valve assembly	Each
 17MP1535-CS	Equinox® II Anti-entrainment ring	Case/10
 17MP1100-90Deg	Equinox® II 90° Scavenger Adapter	Each
 01RE05S10-B	Equinox® II Brass Pin-Indexed Yoke Regulator with 1 x 5/8" Self Sealing DISS Outlets	Each
 17MP1329	Equinox® II 7' replacement hose	Each
 02RT1303	Disposable Mouth Piece (Individually Wrapped)	Case/50
 17MP1997	O-Two 90° (Male to Male) adjustable Oxygen (O ₂) DISS elbow assy	Each
 17MP1998	O-Two 90° (Male to Male) adjustable Nitrous Oxide (N ₂ O) DISS elbow assy	Each
 01CV8028-CS	Single-Use 6 Foot Circuit with Scavenger Hose and Adapter with Mouthpiece	Case/10
 01CV8037-CS	Single-Use Patient Circuit with Scavenger Hose, monitoring line and Adapter with Mouthpiece	Case/10

*Refer to page xx to confirm the hose required for your country

IMMEDIATE CARE



Single Use CPAP

PART N°	DESCRIPTION	UNITS
01CV0212-CS	O-Two Single-Use CPAP System (Adult) c/w Adult #4 Cuff Mask, Head Harness System, CPAP Device with Oxygen Supply Tubing	Case/10
01CV0217-CS	O-Two Single-Use CPAP System (Large Adult) c/w Adult #5 Cuff Mask, Head Harness System, CPAP Device with Oxygen Supply Tubing and Pressure Manometer	Case/10
01CV0218-CS	O-Two Single-Use CPAP System (Adult) c/w Adult #4 Cuff Mask, Head Harness System, CPAP Device with Oxygen Supply Tubing and Pressure Manometer	Case/10
01CV0219-CS	O-Two Single-Use CPAP System (Child) c/w Child #3 Cuff Mask, Head Harness System, CPAP Device with Oxygen Supply Tubing and Pressure Manometer	Case/10
01CV0210-CS	O-Two Replacement Single-Use High Visibility CPAP Flow Vector Valve Only	Case/50



CPAP Accessories

PART N°	DESCRIPTION	UNITS
01CV0220-CS	O-Two Single-Use CPAP Nebulizer System c/w 6 L/min Flow Adapter, TPiece, Updraft Nebulizer with Oxygen Supply Tubing	Case/10
02RT1301-CS	O-Two Single-Use CPAP Nebulizer T-Piece	Case/10
01CV8009-CS	O-Two Replacement Neoprene Head Harness System	Case/5



01BM3201-MO-CS



01BM3211-MO-CS



SMART BAG®

PART N°	DESCRIPTION	UNITS
SMART BAG MO - Disposable		
01BM3201-MO-CS	Smart Bag® w/ manual override control, adult w/ cuffed mask, oxygen reservoir & tubing	Case/12
01BM3211-MO-CS	Smart Bag® w/ manual override control, child w/ cuffed mask, oxygen reservoir & tubing	Case/12
01BM3200-MO-CS	Smart Bag® w/ manual override control, adult w/ Universal Resuscitation mask, oxygen reservoir & tubing	Case/12
01BM3210-MO-CS	Smart Bag® w/ manual override control, child w/ Universal Resuscitation mask, oxygen reservoir & tubing	Case/12
SMART BAG MO - Reusable		
01BM3100-MO	Smart Bag® w/ manual override control, adult w/ silicone mask, oxygen reservoir system in carton	Each
01BM3110-MO	Smart Bag® w/ manual override control, child w/ silicone mask, oxygen reservoir system in carton	Each

SMART BAG® ACCESSORIES

PART N°	DESCRIPTION	UNITS
01BM1000-CS	Adult Ventilation Timing Light, individually packaged with Peel Back Adhesive Backing to attach to any Bag-Valve-Mask Resuscitator	Case/50
01BM1500-CS	Pedi Ventilation Timing Light is individually packaged with Peel Back Adhesive Backing to attach to any Bag-Valve-Mask Resuscitator	Case/50
10-55330-CS	Disposable PEEP Valve 5 - 20 cmH ₂ O	Case/12
01CC8006	Clear Plastic Storage Carrying Case (35mm x 15mm x 14mm)	Each





Mini Ventilation Training Analyzer

PART N°	DESCRIPTION	UNITS
01TA1000	Mini Ventilation Training Analyzer complete with Calibrated Stomach Simulator, Calibrated Lung Simulator, 22/22 mm BVM Connector, in Softpack Carrying Case with Shoulder Strap	Each

Rescue Breather CPR Protective Barrier



PART N°	DESCRIPTION	UNITS
Medical Grade PVC		
02FM4999-CS	O-Two universal mask	Case/12
02FM5000-CS	O-Two Medium Mask (#2) with Anatomical Chin Cup	Case/12
02FM5001-CS	O-Two Medium Mask (#3) with Anatomical Chin Cup	Case/12
Disposable		
02FM5207-CS	O-Two Disposable Cuffed Child/Small Adult #3 Mask with Hook Ring	Case/12
02FM5208-CS	O-Two Disposable Cuffed Large Adult #5 Mask with Hook Ring	Case/12
02FM5210-CS	O-Two Disposable Cuffed Adult #4 Mask with Hook Ring	Case/12






OXYGEN DEMAND VALVES & KITS

Demand Valve and Resuscitator



PART N°	DESCRIPTION	UNITS
01DV2000	O-Two Demand Valve Resuscitator with 6 Foot Supply Hose with 9/16" DISS Fittings & Universal Face Mask	Each
01DV3000	O-Two Demand Valve with 6 Foot Supply Hose with 9/16" DISS Fittings & Universal Face Mask	Each

Demand Valve Accessories

PART N°	DESCRIPTION	UNITS
 01FV4303-DISS	6 Foot Replacement Supply Hose with 9/16" DISS Fittings	Each
 17MP1100-90DEG	Non-Rebreathing Valve Cap Assembly with 90° PEEP Port	Each
 17MP1528	Diaphragm valve assembly	Each
 17MP1535-CS	Anti-entrainment ring	Case/10
 17MP9039	Non-rebreathing valve cap assembly	Each

PRESSURE REGULATORS & OXYGEN ADMINISTRATION

Pressure and Flow Regulators

PART N°	DESCRIPTION	UNITS
Oxygen Regulators - Flow Therapy Only		
01RO01003-B	Therapy Flow 0 - 25 L/min, Barbed Outlet - All Brass Pin-Indexed Yoke	Each
01RO02003-B	Therapy Flow 0 - 25 L/min, Barbed Outlet - All Brass 540 Nut & Stem	Each
Oxygen Regulators - Pressure Only		
01RO01010-B	Brass Pin-Indexed Yoke with 1x9/16" Self Sealing DISS Outlets (No Therapy Flow)	Each
01RO02010-B	Brass 540 Nut & Stem with 1x9/16" Self Sealing DISS Outlets (No Therapy Flow)	Each
01RO02020-B	Brass 540 Nut & Stem with 2x9/16" Self Sealing DISS Outlets (No Therapy Flow)	Each
01RN09N10-B	Brass Pin-indexed Yoke Nitrous Oxide Regulator with 1 X N2O Self Sealing DISS outlet	Each
Oxygen Regulators - Pressure and Flow Therapy		
01RO01023-B	Brass Pin-Indexed Yoke with Therapy Flow 0 - 25 L/min, Barbed Outlet, with 2 x 9/16 Self-Sealing DISS Outlets (Oxygen)	Each
01RO02023-B	Brass 540 Nut & Stem with Therapy Flow 0 - 25 L/min, Barbed Outlet - with 2 x 9/16 Self-Sealing DISS Outlets (Oxygen)	Each





Cylinders

PART N°	DESCRIPTION	UNITS
01EQ7007	Entonox® "D" Cylinder c/w Pin Indexed Medical Post Valve	Each
01EQ7008	Entonox® "E" Cylinder c/w Pin Indexed Medical Post Valve	Each

Pressure Regulators & Oxygen Administration Accessories

PART N°	DESCRIPTION	UNITS
01OR1002-CS	Viton/Brass Bodoc Post Valve Gasket (Yoke Seal)	Case/10
01OR1007	Low Profile High Pressure Contents Gauge (0 - 3,000 psi)	Each
02RT5400-10	Nylon Oxygen Cylinder Wrench with PVC Strap	Case/10
01OR6664-CS	Oxygen Nut & Nipple Connector	Case/10




MOUNTING AND CARRYING SOLUTION

Equipment Carrying Cases

PART N°	DESCRIPTION	UNITS
 01CV7050	eSeries Automatic Transport Ventilator Carrying Case - With sling-style shoulder strap (specifically designed for eSeries)	Each
 01CC8005	Emergency Responder Moulded Plastic Carrying Case in International Orange for "D" or "Jumbo D" Cylinders	Each
 01CC8017	Softpack Clamshell Carrying Case with Sling Style Shoulder Strap for "D" or "Jumbo D" Cylinders	Each

Ventilator Mounting and Carrying Systems

PART N°	DESCRIPTION	UNITS
 01CV7035	O-Two Medical "e" Series Ventilator "Smart Mount" Multi-configuration Mounting Bracket (Ambulance Cot, Hospital Stretcher, Bed, Roll Stand)	Each



www.otwo.com

45A Armthorpe Road, Brampton, ON, Canada, L6T 5M4

Telephone: +1 905 792-OTWO (6896)

N.A. Toll Free: +1 800 387 3405 Facsimile: +1 905 799 1339

Email: resuscitation@otwo.com

o_two controlled[™]
ventilation