

Oxygen *Flow* Controllers



The O-Two adjustable oxygen therapy flow controller offers reliable performance with high-quality manufacturing. It's a 12-position, calibrated, fixed-orifice flow controller that operates in any position. Available with 9/16" female and male DISS inlets and 9/16" male outlets.

Lightweight and ergonomically designed, it's compatible with humidifiers, making it ideal for adjustable oxygen therapy applications with superior quality, size, and performance.

Oxygen Therapy Treatment for Single Patient Incidents

SPECIFICATIONS

THERAPY FLOW OPTION 1	0, 0.5, 1, 2, 3, 4, 5, 6, 8, 10, 12, 15 (+/- 10%)
THERAPY FLOW OPTION 2	0, 1, 2, 3, 4, 6, 8, 10, 12, 15, 20, 25 (+/- 10%)
MAXIMUM RATED INPUT PRESSURE	50 PSI
OPERATING TEMPERATURE RANGE	-18°C to 50°C (0°F to 122°F)
STORAGE TEMPERATURE RANGE	-40°C to 60°C (-40°F to 140°F)
SIZE (EXCLUDING OUTLET AND INLET FITTING) (H X W)	2.0 x 1.37 in / 5.1 x 3.5 mm
WEIGHT	4.0 oz / 0.1 kg

ORDERING INFORMATION

01FC0815D

Adjustable 0 - 15 L/min Click Stop Flow Controller - 9/16" DISS Oxygen input female connector (Available in Canada Only)

Each

01FC0825D

Adjustable 0 - 25 L/min Click Stop Flow Controller - 9/16" DISS Oxygen input female connector (Available in Canada only)

Each

01FC0115D

Adjustable 0 - 15 L/Min Click Stop Flow Controller - 9/16" DISS Oxygen input male connector

Each

01FC0125D

Adjustable 0 - 25 L/Min Click Stop Flow Controller - 9/16" DISS Oxygen input male connector

Each

01OR6664-CS

Oxygen Nut & Nipple Connector

Case/10

www.otwo.com

45A Armthorpe Road, Brampton, ON, Canada, L6T 5M4
 Telephone: +1 905 792-OTWO (6896) N.A. Toll Free: +1 800 387 3405
 Facsimile: +1 905 799 1339 Email: resuscitation@otwo.com



MedNet EC-REP llb GmbH Borkstrasse, 10
 48163 Münster, Germany

