

# Disposable Circuit with Mask



**O-TWO EQUINOX® RELIEVE  
ANALGESIC GAS DELIVERY SYSTEMS  
SINGLE-USE PATIENT CIRCUITS**  
**HIGH QUALITY/PERFORMANCE PATIENT  
CIRCUITS FOR GAS DELIVERY WITH MASK**

Designed for use with any of the currently available gas delivery systems, these high-quality, high-performance single patient-use circuits are an essential adjunct to providing low work of breathing and ease of patient management.

The single patient-use circuits incorporate both inspiratory and scavenging port with a expiratory one-way valve to ensure fresh gas delivery and no re-breathing of exhaled gas. The scavenging port is compatible with any universal 30 mm scavenger hose available in the market, eliminating the need to purchase special hoses.



## SPECIFICATIONS

<b>PATIENT VALVE (L X W X H):</b>	
<b>INCHES:</b>	4.5 x 2.5 x 2.0 in
<b>CENTIMETRES:</b>	11.4 x 6.3 x 5 cm
<b>TUBING LENGTH COLLAPSED:</b>	12 in / 30.4 cm
<b>TUBING LENGTH EXPANDED:</b>	36 in / 91.4 cm
<b>WEIGHT</b>	2.8 oz (80 g)
<b>OPERATING TEMPERATURE RANGE:</b>	-18 to 50°C (0 to 122°F)
<b>STORAGE TEMPERATURE RANGE:</b>	-40 to 60°C (-40 to 140°F)
<b>RH:</b>	15 - 95%

## ORDERING INFORMATION

**01CV8019** - O-Two Gas Delivery Systems Single Use Patient Circuits with Universal Mask



[www.otwo.com](http://www.otwo.com)

45A Armthorpe Road, Brampton, ON, Canada, L6T 5M4  
Telephone: +1 905 792-OTWO (6896) N.A. Toll Free: +1 800 387 3405  
Facsimile: +1 905 799 1339 Email: [resuscitation@otwo.com](mailto:resuscitation@otwo.com)